



FT 13 Thermometer



beurer 37.0°C



Flexible measuring tip



Audible warning signal if temperature exceeds 37,8°C (fever alarm)



Waterproof

Digital thermometer

Measurements in approx. 30 seconds

Without glass, mercury free

1 memory

Measurement accuracy:

±0.1°C: 35.5°C to 42.0°C (±0.2°F: 95.9°F to 107.6°F)

±0.2°C: 32.0°C to 35.4°C; 42.1°C to 42.9°C (±0.4°F: 90°F to 95.7°F; 107.8°F to 109,9°F)

Audible signal when measuring stops

Extra large display

Can be disinfected

With storage cap

Incl. 1 x 1,5V Batterie LR41

Battery changing indicator

Automatic switch off

Delivered in sales display

Indication of °C and °F

Product measurements (LxWxH): 142 x 23 x 13 mm

Product weight: approx. 12,9 g (incl. Batterie)

Medical device

5 year guarantee

Sales unit: 20 / Shipping carton: 10 x 20

EAN no.: 4211125 79109 4

Item no.: 791.09

791.09_FT13_2018-04-18_01_DS_EN_BEU Subject to errors and changes