



HK 145
Cosy Weight
beurer

HK 145 Cosy Weight Green Planet




Upper side / underside from 100 % recycled material



Control housing with approx.
50 % recycled material



Packaging of 100 %
recycled material



Self-hugging warmth thanks to the wonderfully snuggly heat pad

Can be used on different areas of the body, e.g. stomach, back, neck, shoulders

Optimised feeling of warmth

Long-lasting heat release thanks to glass bead filling

6 illuminated temperature settings

Breathable, cosy and skin-friendly

Electronic temperature control

Fast heating

Automatic switch off after approx. 90 minutes

BSS (Beurer Safety System)

Detachable control unit

Machine washable up to 30 °C extra gentle cycle

OEKO-TEX®

The textiles used for this device meet the stringent human ecological requirements of OEKO-TEX® Standard 100, as verified by Hohenstein Research Institute.

Mains operation:

220-240 V; 50-60 Hz; 100 W

Product dimensions (LxB): 60 x 30 cm

Product weight: approx. 1500 g

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 421112520039 8

Item no.: 200.39



Top and bottom contain 100 %
GRS certified recycled polyester
Control Union 1142381

200.39_HK145-GP_2024-08-16_02_DS-EN-BEU Subject to errors and changes