




BF 183 Diagnostic bathroom scale

 **german**
engineering


With extra large display



Diagnostic scale with visceral fat and metabolic age

Digit size 38 mm

Weighing area made from safety glass (6 mm)

Electrodes made of brushed stainless steel

Determination of body weight, fat, water, muscle density, BMI, visceral fat and metabolic age

Capacity 180 kg

Introduction:

- 0.1% for body fat, body water and muscle density

- 100 g for body weight

Change over to kg/lb/st

10 user memories

Automatic switch-on and switch-off

Overload indication

1 x 3V CR 2032 battery included

Product measurements: 30 x 30 x 1,8 cm

Product weight: approx. 1549 g

5 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 10054 4

Item no.: 100.54