



GL 44  
mmol/L

**beurer**  
medical

## GL 44 mmol/L Blood glucose monitor

 **german  
engineering**



With USB cable and Health-Manager software download for free



Illuminated insertion slot



White illuminated XXL display



Starter Pack: measuring unit, lancing device, 10 lancets, 10 test strips, batteries, pouch, instruction manual, short instruction, blood glucose diary, USB cable

Only 0.6 µl of blood necessary

Period of measurement: approx. 5 seconds

480 stored results with date, time (data memory also with battery change)

Stored results can be scrolled forward and back

Average value calculation for 7, 14, 30 or 90 days for  
- all measurements                      - marked values

Blood sample: fresh whole blood (capillary)

Haematocrit: 20 - 60%

Plasma calibrated

Measurement range: 1,1 - 35,0 mmol/L

Blood volume check: measurement only possible with sufficient sample of blood

Complies with EN ISO 15197:2015 in its entirety (obligatory since Mai 2016)

Automatic switch off

Product measurements: 52 x 95 x 16 mm

Product weight: approx. 44 g (including batteries)

2 x Lithium batteries CR 2032 for more than 1.000 measurements

Storage/Transport:	Temp.:	Rel.humidity:
Unit:	-10°C to +65°C	< 90 %
Test stripes:	+2°C to +30°C	< 90 %
Operating:	+10°C to +40°C	< 90 %

Medical product

5 year guarantee

	Black	Purple	White
Item no.:	463.16	463.14	463.12
EAN no.:	4211125 463169	4211125 463145	4211125 463121

GL44\_mmol\_L-2019-01-31-01\_DS\_EN Subject to errors and changes



GL 44  
mmol/L

**beurer**  
medical



Patented\* soft & safe lancing device




Illuminated test strip insertion slot and wide test strips



Simple measurement marking for optimised diabetes management: before eating, after eating, general marking



## CONNECTION



GL44\_mmoL\_2019-01-31\_01\_DS\_EN Subject to errors and changes

\* tested online application "beurer HealthManager" Cloud