



BY 80 Baby scale


**german
engineering**

**LCD
display**

858

 Easy-to-read LCD display
72 x 28 mm

**20 kg
capacity**

20 kg capacity


**hold
function**

 Automatic and manual hold
function for recording the weight
of wriggling babies


Automatic and manual hold function for recording the weight of wriggling babies

Large LCD display (digit size 23 mm)

Curved weighing surface

Switch-on technology: button-on

Weight capacity: 20 kg

Graduation: 5 g

Adjustable between kg/lb/oz

Easy to clean

Non-slip, abrasion-resistant rubber feet

Automatic switch-off

2 x 1.5V AA batteries

5-year guarantee

Pack: 2/Shipping carton: -

EAN no.: 4211125956059

Item no.: 956.05