



## IH 58 Kids Nebuliser




compressed  
air technology

Nebulisation of liquid medicine using compressed air technology



focused on  
babies and kids

Effective inhalation from the 6th month of life\*



nebulisation  
performance

Nebulisation performance: approx. 0.25 ml/min



For the treatment of the upper and lower respiratory tract

For use in colds, asthma, respiratory disease, etc.

- for fast and effective inhalation

High proportion of respirable particles

High nebulisation capacity (approx. 0.25 ml/min.)

- short inhalation time

Extensive accessories:

With a removable "giraffe", mouthpiece, silicone children's and baby mask, angled fitting, atomiser, 10 filters, hose and storage bag

Can be disinfected

Operating pressure: approx. 0.25-0.5 bar

Medical device

Mains operation:

Input: 100-240 V~; 50-60 Hz; 0.5 A

Output: 5 V === 2 A

Product measurements (LxWxH): 110 x 62 x 47 mm

Product weight: approx. 176 g

5 year guarantee

Sales unit: 3 / Shipping carton: 4 x 3

EAN no.: 4211125 60217 8

Item no.: 602.17



\*Inhalation depends on the fit of the mask.